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 **NOTE:** For U. S. customers please use language column marked with "US".

1 Intended use

Light cables serve to transmit the light from the external light source to the site of the intervention during diagnostic and therapeutic interventions in endoscopy and exoscopy.

Contraindications: No contraindications directly related to the product are currently known.

2 Safety information

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of infection: These light cables are not sterile when delivered. The use of non-sterile light cables poses a risk of infection for patients, users and third parties. Inspect the light cables for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the lights cables before initial use and before and after every subsequent use using validated procedures.

WARNING: Risk of injury: Incorrectly assembled and damaged medical devices can lead to injuries to the patient. The medical devices and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

3 Handling

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

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
WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope. Always wear standardized protective equipment when using lasers.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

- Avoid direct tissue contact with the distal end, as well as with the light connection of the telescope and the light cable.
- Do not rest the endoscope or the light cable on the patient or where they are in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a type CF applied part is used. In this case, an endotherapy device with a type CF applied part must be used in order to reduce the total patient leakage current.

CAUTION: Damage to the light cable! Light cables must not be kinked or squeezed. A minimum radius of 10 cm must be ensured when coiling the cable as the sheath or the optical fibres may otherwise be damaged.

- Avoid squeezing, kinking, stretching, or sharply bending the cable.
- Secure the light cable to the telescope with a 1/4 turn (adaptors are available to connect KARL STORZ light cables to other manufacturer's telescopes and light sources).
- To mount the light cable **495 NVC**, first remove the KARL STORZ light cable adaptor 495 G from the light connection. Then attach the light cable directly to the telescope thread (Fig. 1, 2).
- Hold the light cable by the connector handle when connecting it to or disconnecting it from the light source
- Avoid puncturing the external sheath with towel clips. Never squeeze the cable between forceps.

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* If there is an identification symbol on the light cable and the endoscope (see table right), assign matching combinations according to the symbols shown.

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3.1 Overview: Recommended combinations*

Light cable diameter		Endoscope diameter	
4.8 – 5.0 mm		6.6 – 12.0 mm	
3.0 – 3.5 mm		3.0 – 6.5 mm	
2.0 – 2.5 mm		0.8 – 2.9 mm	

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4 Reprocessing

! WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

! WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

! WARNING: If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.

! CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiologic range of action of the chemicals used.

! CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the Internet at www.karlstorz.com.

! CAUTION: National laws and regulations must be observed.

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i **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

4.1 Accessories

Accessories required for carrying out reprocessing:
Brushes: 27652

4.2 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, pre-clean the medical device, e.g., by wiping it down and rinsing it. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

4.3 Manual cleaning

The medical device must be completely immersed in a cleaning solution. After the necessary exposure time, clean the instrument mechanically with the aid of brushes or a sponge. Finally, it must be rinsed with cold water to ensure neutralization.

4.4 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints and openings are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

4.5 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Automated cleaning/chemical-thermal disinfection (only applies to fluid light cables)

A maximum temperature of 65 °C must not be exceeded. The relevant national regulations and the chemical manufacturer's instructions must be taken into account when using this procedure.

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**Machine cleaning/thermal disinfection
(does not apply to fluid light cables)**

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off afterwards by hand.

4.6 Inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Threads and glide surfaces must be lubricated locally with instrument oil.
- Afterwards, a functional check must be carried out.

4.7 Inspecting the light cables

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection. Such dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the light intensity. As of a certain number of broken optical fibers, the light intensity deteriorates considerably. Depending on the application and the combination of light source, light cable and applied part used, the light intensity is then no longer sufficient. Check that the available light intensity still ensures safe working prior to use. If this is not the case, the fiber optic light cable should no longer be used.

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4.8 Drying the optical end faces

The fiber areas in the light inlet piece in particular must be dried carefully with 70 % alcohol. Disinfectant and cleaning agent residue, e.g., in the light inlet piece, can become burnt on when the light cable is connected, thereby seriously impairing light transmission.

- NOTE:** The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).
- NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques' (Art. no. 96211004).

4.9 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

4.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

4.10.1 Steam sterilization using the fractionated prevacuum procedure (does not apply to fluid light cables)

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132 °C–137 °C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

- NOTE:** Greased components must be disassembled for sterilization, so that steam penetration is ensured.
- NOTE:** This validation is applicable for reusable tube sets with a maximum length of 200 cm.

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4.10.2 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30–45 minutes.

WARNING: Medical devices must be dismantled for sterilization.

CAUTION: For gas sterilization with ethylene oxide, the airing times for instruments stipulated by the device manufacturer must be observed due to gas absorption of the materials. The airing times depend on the processing technology of the EO sterilizer (concentration, process control).

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

4.10.3 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Medical devices must be dismantled for sterilization.

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

CAUTION: For detailed information, please consult the user handbook of the respective device.

NOTE: Please refer to the „STERRAD® Sterility Guide“ to verify whether the respective medical device can be sterilized in the various STERRAD® devices.

NOTE: Sterilization is not possible on surfaces which have been greased and lubricated.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 50, 100S, 200 without booster
- STERRAD® NX® Standard Cycle
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Standard Cycle
- STERRAD® 100NX® Flex Cycle

4.10.4 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™ 1

Detailed information on the selection of the appropriate cycle for the various device generations is available from the manufacturer of STERIS®.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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4.11 Limits of reprocessing

The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

4.11.1 Fluid light cable with sterilization method STERRAD® 100NX®

Provided the application and reprocessing are done appropriately and professionally, KARL STORZ validated 50 cycles.

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5 Operation and storage conditions

	Temperature	Rel. humidity
Storage	-20°C ... +60°C	10% ... 90%
Application	+10°C ... +40°C	30% ... 70%
Air pressure	700 hPa – 1060 hPa	

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6 Directive compliance

This medical device bears a CE mark in accordance with MDD 93/42/EEC (does not apply to devices with article numbers 69495 xx).

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7 Disposal

For disposal, no special measures are necessary. Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

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8 Packaging symbols

LOT	Batch code
REF	Catalogue number
QTY	Number of products in the product packaging
	Fragile, handle with care
	Temperature limit

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1 Intended use

Light cables are used to conduct light from a light source to an endoscope.

Contraindications: No contraindications relating directly to the product are currently known.

2 Safety information

WARNING: Risk of injury and damage to the products: Disregarding these user instructions and all of the user instructions for the products used in combination with one another can lead to injuries in patients, users and third parties or to damage to the products themselves. Please read all relevant instruction manuals thoroughly and always follow all instructions. Check the functionality of the products used in combination.

WARNING: Risk of infection: These instruments are not delivered in a sterilized condition. Using non-sterilized instruments bears a risk of infection for patients, users and third parties. Check instruments for visible contaminations. Visible contaminations indicate missing or incorrect preparation. Prepare the instruments prior to first use as well as before and after each subsequent use; apply validated procedures.

3 Handling

WARNING: Risk of injury: Improper use of medical equipment bears a risk of infection for patients. Users of medical equipment have to be sufficiently medically qualified and familiar with the use of the respective equipment.

WARNING: Risk of injury: A capacity overload due to excessive application of force can result in cracks, bending or malfunctions of the medical equipment and lead to injuries in the patient or user. Do not attempt to bend back twisted equipment into its original shape.

WARNING: Never place the end of a fiber optic light cable or telescope connected to a light cable on or under a surgical drape while the light source is on. The intensity of the light source may cause burns to the patient and/or the surgical drape. KARL STORZ recommends that the appropriate size light cable be used to minimize the heat energy generated from the light source. Use of appropriate size light cable will reduce the risk of patient burns, as well as the risk of igniting flammable material. Never attach a cable or cable connector to patient or drape.

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WARNING: Risk of injury: Optical radiation can lead to eye injuries. Never look into the free end of a connected light conductor or endoscope. Always wear proper standard protective gear during use of lasers.


WARNING: Risk of injury and damage to the products: The high light energy could heat up the distal end, the light connection, as well as adjacent components and material in front of the light emission window. This can lead to burns in the patient, the user or the surgical equipment.

- Avoid direct material contact with the distal end and the light connection of the endoscope and the light cable.
- Do not leave the endoscope and/or light cable on the patient or in direct contact with surgical equipment.
- Always select a light output for the cold light source that is as low as possible to barely achieve optimal illumination of the surgical field.

WARNING: Risk of injury: Patient leakage currents may add up when energy-supplied endoscopes and endotherapy devices are used simultaneously. This is especially important if an endoscope is used with an application part of type CF. In such a case an endotherapy device with an application part of type CF has to be used in order to minimize the total patient leakage current.

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- Avoid squeezing, kinking, stretching, or sharply bending the cable.
- Secure the light cable to the telescope with 1/4 turn (Adapters are available to connect KARL STORZ light cables to other manufacturer's telescopes and light sources).
- To connect the **495 NVC light cable**, first remove the KARL STORZ light cable adaptor 495 G from the light connection. Then attach the light cable directly to the telescope thread (Fig. 1, 2).
- Grasp the cable by the strain relief handle when inserting the cable into or removing it from the light source.
- Avoid puncturing the external sheath with towel clips. Never squeeze the cable between forceps.

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Diameter light cable		Diameter endoscope	
4.8mm – 5.0mm		6.6mm – 12mm	
3mm – 3.5mm		3mm – 6.5mm	
2mm – 2.5mm		0.8mm – 2.9mm	

4 Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing.

CAUTION: National laws and regulations must be observed.

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4.2 Pre-Cleaning

1. Immediately after use, at the point-of-use (e.g. operating room), wipe the ligh cable with a soft, lint-free disposable cloth moistened with water or a diluted mild/neutral pH enzymatic cleaner (e.g. Enzol) to remove all heavy soiling.

2. Contain and transport the contaminated light cable to the decontamination area for cleaning.

NOTE: If immediate cleaning is not feasible, place the endoscope in a container and soak for no longer than 30 minutes in diluted neutral/mild pH enzymatic cleaning solution (e.g. Enzol).

4.3 Manual cleaning

1. Prepare the cleaning solution per manufacturer's instructions.
2. Thoroughly rinse the endoscope and light cable for a minimum of 2 minutes with water to remove all gross debris.
3. Using a low-lint soft cloth dampened with cleaning solution, wrap the cloth around the light cable and wipe the entire cable portion from end to end, back and forth, for a minimum of two passes.
4. Completely immerse the light cable in the cleaning solution and keep immersed for a minimum of 5 minutes.
5. While immersed, brush the hard surfaces of the light cable with a soft bristle brush. Pay close attention to the stainless steel locking mechanism, turning the locking mechanism left and right, brushing for a minimum of 30 seconds. Do not brush the fiber optic ends of the light cable. Instead, use a lint-free, soft cloth to wipe the fiber optic ends of the cable and all flexible areas of the cable. Rinse the bristles of the brush with water if heavy debris collects on it. Keep brushing until there is no visible debris on the bristles of the brush (minimum of 2 minutes).
6. After the immersion period, remove the light cable from the cleaning solution and rinse with cold running tap water for a minimum of one minute.
7. Wipe the glass window of the light cable with a dampened cloth.
8. Dry the light cable with a soft, lint-free cloth or clean, dry, oil-free compressed air.

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4.4 Inspection and Care

The cleaned medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning once more.
- Damaged or corroded medical devices must be withdrawn from use.

4.5 Inspecting the Light Cables

Hold one end of the light cable in the direction of a light source and check the number of dark points at the other end. These dark points represent broken optical fibers in the light guide bundle. Individual broken optical fibers do not necessarily mean a noticeable reduction in image quality. Depending on the combination of light source, light cable and application part, light transmission is no longer sufficient. The light cable should no longer be used then.

4.6 Packaging Systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

4.7 Sterilization

KARL STORZ has validated the following sterilization methods for the light cables:

1. Fiber Optic Light Cable:

- Steam Sterilization
- Ethylene Oxide Sterilization
- STERRAD® Sterilization Systems
- AMSCO® V-PRO™ Sterilization Systems

2. Fluid Light Cable

- STERRAD® NX Sterilization System
- STERIS® System 1E Liquid Chemical Processing System

Refer to the appropriate section for full instructions for use and parameters for each sterilization methods.

4.8 Steam Sterilization

The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following steam sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

1. Place the light cable in a steam sterilization tray or container.
2. Double wrap the sterilization tray and place in the sterilizer.

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3. Use the following parameters for sterilization:

Temperature:	270°F (132 °C)
Exposure Time:	4 minutes
Dry Time:	20 minutes

4.8.1 Ethylene Oxide Sterilization (EO)

1. Place the light cable in a EO gas sterilization tray or container.
2. Double wrap the sterilization tray and place in the sterilizer.
3. Process in 100% Ethylene Oxide Sterilizer under the following parameters:

Conditioning Parameters

Temperature:	55 °C
Humidity:	70% RH
Conditioning Dwell Time:	30 minutes

Sterilization Parameters

Temperature:	55 °C
Humidity:	70%
Humidity Dwell Time:	30 – 45 minutes
EO Gas Concentration:	725 ± 30 mg/L
EO Gas Exposure Time:	60 minutes

Aeration Parameters

Time:	12 hours
Aeration Temperature:	55 – 59 °C

**4.8.2 Hydrogen peroxide (H₂O₂)
 sterilization – ASP STERRAD®**

WARNING: Medical devices must be disassembled for sterilization.

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect to fluid light cables.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for these medical products:

- STERRAD® 100S
- STERRAD® NX® Standard Cycle
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Standard Cycle
- STERRAD® 100NX® Flex Cycle

CAUTION: Fluid Light Cables may only be processed in the STERRAD® NX® Sterilization System using the Advanced Cycle.

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1. Place the light cable in a STERRAD® compatible sterilization tray or container.
2. Double wrap the tray and place in the sterilizer.
3. Start the sterilizer using the appropriate cycle.

**4.8.3 Hydrogen peroxide (H₂O₂)
 sterilization – STERIS® AMSCO®
 V-PRO™ 1**

The following AMSCO® V-PRO™ sterilization processes have been validated and approved by KARL STORZ for these medical products:

- V-PRO™ 1 Plus Lumen and Non-Lumen Cycle
- V-PRO™ maX Lumen and Non-Lumen Cycle

1. Place the light cable in a V-PRO™ compatible sterilization tray or container.
2. Double wrap the tray and place in the sterilizer.
3. Start the sterilizer using the appropriate cycle.

**4.8.4 Chemical Low-Temperature
 Sterilization with Peracetic Acid –
 STERIS® SYSTEM 1E**

1. Place the light cable in the appropriate processing tray.
2. Load the tray into the SS1E and start the system.

 **NOTE:** Refer to the SS1E operator's manual for complete instructions for use.

4.8.5 High Level Disinfection

 **WARNING:** High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

 **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time, and service life.

This device may be chemically disinfected using CIDEX Activated Dialdehyde, a 2.4% concentration of glutaraldehyde solution, or CIDEX OPA, a 0.55% concentration of ortho-Phthalaldehyde solution.

1. Prepare the disinfecting solution per manufacturer's instructions.
2. Place the device into the disinfection solution. Ensure that no air bubbles are present on the surface of the device.
3. Use the following KARL STORZ validated conditions to achieve high-level disinfection:

CIDEX Activated Dialdehyde: Immerse for 45 minutes at 77°F (25°C).

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- After disinfection is complete, remove the device from the disinfection solution and completely immersed in water. Keep immersed for 1 minute.
- Discard the water and repeat rinse for an additional 2 immersion rinses.
- Dry the light cable with a soft, lint-free, sterile cloth or clean, dry, oil-free compressed air.

4.9 Reprocessing Limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

5 Application/Storage Conditions

	Temperature	Relative humidity
Storage	-4°F (-20 °C) ... 140°F (60 °C)	10 % ... 90 %
Application	50°F (10 °C) ... 104°F (40 °C)	30 % ... 70 %
Air pressure	700 hPa–1060 hPa	

6 Directive Compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

7 Disposal

To dispose of fiber optic light cables, no special measures are necessary. National laws and regulations must be observed.

8 Packaging symbols

	Batch code
	Catalogue number
	Number of products in the product packaging
	Fragile, handle with care
	Temperature limit
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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