

Aura KTP

Laser Fiber



Medical Device ACCESSORIES

The Aura KTP Laser Fiber is a buffered 3.5 meter long fiber with a flat tip and stripped end. The laser fiber may be used for a variety of surgical procedures and is only intended for use with an Aura KTP laser.

- Integrated adapter for ease of connection with the Aura KTP Laser
- Stripped and pre-cleaved for convenience
- Packaged sterile for single-use
- Available in 400 (0.4mm) and 600 (0.6mm) micron sizes
- Stripper and cleaver tools available with clear sterilization instructions



Core Size	BSC Part #	ForTec Part #
0.4mm/400µm	10-0622	KTP-400
0.6mm/600µm	10-0612	KTP-600



[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

to start ForTec Mobile Service