



Case-by-Case
MOBILE SERVICE

Aura XP™

Laser System

The Aura XP laser system delivers steady-state reliability that can be applied to a wide range of soft tissue clinical applications such as delicate laryngeal situations.

- Precise Cutting & Coagulation
- Unique Starpulse™ Technology
- KTP pulsed 532nm wavelength
- Selection of handpieces for varied procedures and situations
- Range of fiber sizes for diversity and surgeon preference



[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

to start ForTec Mobile Service

AuraXP™ Laser System

Laser Specifications:

Laser type: Quasi-continuous KTP 532nm

Power: 0.5 to 15W

Dimensions: 16" (h) x 12" (w) x 23" (d)

Weight: 60 lbs.

Precise Cutting & Coagulation

Starpulse™ technology precisely limits the depth of penetration for increased control, while providing excellent hemostasis, vaporization and cutting power.

Unique Starpulse™ Technology

Delivers a pulse width as short as 1 millisecond to minimize thermal spread.

KTP pulsed 532nm wavelength

Highly absorbed by oxyhemoglobin, delivering precise coagulation while leaving the surrounding tissue unaffected.

Selection of handpieces for varied procedures & situations

ENT/Laryngeal – straight and curved [as shown]

Sinus – straight, curved right and curved left

Range of fiber sizes for diversity & surgeon preference

200, 300, 400 and 600 micron

