



Case-by-Case
MOBILE SERVICE

EVOLVE[®]

180W High Powered Diode

The EVOLVE[®] Laser is used to vaporize and remove the obstructing prostate tissue that causes the symptoms associated with BPH. The 980nm wavelength provides the highest combined absorption in water and hemoglobin. Absorption of the laser energy by hemoglobin leads to excellent hemostasis. This 980nm laser allows for precise ablation of the obstructing tissue making it an ideal wavelength for this procedure.



Laser Specifications:

Laser type: Diode laser
Wavelength: 980nm
Power: up to 180W
Dimensions: 12" (h) x 22" (w) x 12" (d)
Weight: 66 lbs.

[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

www.fortecmedical.com