



FORTECmedical

>> *FORWARDING* technology



Case-by-Case MOBILE SERVICE

ShockPulse-SE

Lithotripsy System

The ShockPulse-SE (Stone Eliminator) is the next generation Cyberwand system from Olympus, offering dual action lithotripsy with proprietary technology to simultaneously produce constant ultrasonic wave energy and intermittent shockwave energy via a revolutionary single probe design. The ShockPulse offers unmatched performance and greatly reduces procedure time.



- **Superior Speed & Efficiency:**
Effectively fragments and aspirates all stone sizes, shapes and composition
- **Intuitive Surgeon Control:**
Ergonomic hand piece with integrated suction and hand activation
- **Versatile:**
Full range of probes and sizes (kidney, ureter & bladder)

[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

to start ForTec Mobile Service

ShockPulse-SE


Lithotripsy System

Laser Specifications:

Voltage: 90 - 264 VAC
Frequency: 50/60 Hz
Dimensions: 14.25" (w) x 15.5" (h) x 4.25" (d)
Weight: 14.1 lbs.
Transducer Dimensions: 7.7" (l) x 2" (dia.)
Cable Length: 10 ft.


Intuitive Surgeon Control Hand Activation and Suction

- Single hand piece design, no need to connect pneumatic and ultrasonic transducers
- New transducer hand activation eliminates the need for a footswitch
- Ergonomic design allows for streamlined control during PCNL, URS and Bladder procedures
- Transducer maintains lower surface temperature throughout the entire procedure



Probe Options

Procedure	Probe Size	Hub Color
PCNL	3.76mm	Silver
	3.4mm	Blue
Mini PCNL, PCNL	1.83mm	Red
URS	1.5mm	Green
	.97mm (no suction)	Gold
Bladder	3.76mm	Silver
	3.4mm	Blue



FORTECmedical

>> Start **FORTEC** Mobile Service **Today**

www.fortecmedical.com | 800.963.7101