



**Case-by-Case
MOBILE SERVICE**

Single Use
Morcellator

The Cyber Blade™ single use morcellator provides a convenient all-in-one solution for individual procedures. The kit includes a complete sterile morcellator device: blade and handpiece (with motor compartment integrated), suction tube and tissue collector.

Cordless Convenience

- All functions are hand controlled
- Handpiece is light and ergonomic
- No need for motor and vacuum console



[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

to start ForTec Mobile Service

Single Use Morcellator

Specifications:

Classification: Battery-powered system (rechargeable)

Morcellator Dimensions: 24" x 2.2" x 1.5"

Blade Length: 14.5"

Blade Diameter: 13.5Fr

Power Supply: 110/220 VAC; 50/60Hz

Advantages of a Single Use Instrument



Eliminate reprocessing issues



A sharp new high-performance blade for every case



No cross-contamination risk among cases



Eliminate maintenance and repair costs



No delays due to instrument issues



Fixed cost per treatment

