



FORTEC_{medical}

AuraXP™

Laser System



Case-by-Case
MOBILE SERVICE

The Aura XP laser system delivers steady-state reliability that can be applied to a wide range of soft tissue clinical applications such as delicate laryngeal situations.

Precise Cutting & Coagulation

Unique Starpulse™ Technology

KTP pulsed 532nm wavelength

Selection of handpieces for varied procedures and situations

Range of fiber sizes for diversity and surgeon preference



✓ **Advanced Surgical Technologies.**

✓ **Clinical Procedural Support.**

✓ **Access for All.**

Product Features & Specifications

AuraXP™ Laser System

Precise Cutting & Coagulation

Starpulse™ technology precisely limits the depth of penetration for increased control, while providing excellent hemostasis, vaporization and cutting power.

Unique Starpulse™ Technology

Delivers a pulse width as short as 1 millisecond to minimize thermal spread.

KTP pulsed 532nm wavelength

Highly absorbed by oxyhemoglobin, delivering precise coagulation while leaving the surrounding tissue unaffected.

Selection of handpieces for varied procedures & situations

ENT/Laryngeal – straight and curved [as shown]
Sinus – straight, curved right and curved left

Range of fiber sizes for diversity & surgeon preference

Range of fiber sizes for diversity & surgeon preference

Laser Specifications:

Laser type: Quasi-continuous KTP 532nm

Power: 0.5 to 15W

Dimensions: 16" (h) x 12" (w) x 23" (d)

Weight: 60 lbs.

