



FORTECmedical

# CO<sub>2</sub>RE Fractional CO<sub>2</sub> Resurfacing System



Case-by-Case  
MOBILE SERVICE

The fractional CO<sub>2</sub> laser system is capable of targeting multiple dermal levels as well as performing traditional skin resurfacing, ablation and excision of lesions. The novel technology also provides a distinct set of pattern generation abilities, increasing the range of possible treatments.

**Unrivalled Versatility**

**Exceptional Control**



✓ **Advanced Surgical  
Technologies.**

✓ **Clinical Procedural  
Support.**

✓ **Access for All.**

## Product Features & Specifications

# CO<sub>2</sub>RE Fractional CO<sub>2</sub> Resurfacing System

### Unrivaled Versatility

Seamless transitions between a variety of different modes enable more treatment options for a wider scope of patient needs.

This laser can also be turned for superficial mid-level, or deep dermal procedures. CO<sub>2</sub>RE Fusion technology even allows for the simultaneous treatments of multiple levels.

#### Laser Specifications:

**Laser type:** CO<sub>2</sub>

**Wavelength:** 10.6  $\mu$ m

**Power:** Up to 60W

**Dimensions:** 70" (h) x 17.3" (w) x 15.75" (d)

### Exceptional Control

An intuitive touch-screen displays mode settings and enables simple switches between and control over intensity, power and depth. In addition, a thin, lightweight handpiece is ideal for optimal maneuverability and precision during procedures.

### Reusable Handpiece

