

Morcellator

Single Use



Case-by-Case
MOBILE SERVICE

The Cyber Blade™ single use morcellator provides a convenient all-in-one solution for individual procedures. The package includes a complete sterile morcellator device: blade and handpiece (with motor compartment integrated), suction tube and tissue collector.

Cordless Convenience

All functions are hand controlled

Handpiece is light and ergonomic

No need for motor and vacuum console



Optimized Design

- oscillating blade movement
- window designed to minimize tissue detachment

Lever for Suction Adjustment

Blade Activation Buttons

Device Activation & Adjustment of Blade Speed

Reusable Internal Unit

- high-performance motor
- rechargeable battery

Vacuum Coupling

- No strict need for an external vacuum station: if not available, connect the handpiece directly to O.R.'s vacuum splitter.

Product Features & Specifications

Morcellator Single Use

Advantages of a Single Use Instrument



Eliminate reprocessing issues



No cross-contamination risk among cases



No delays due to instrument issues



Eliminate maintenance and repair costs



A sharp new high-performance blade for every case



Fixed cost per treatment

Laser Specifications:

Classification: Battery-powered system (rechargeable)

Morcellator Dimensions: 24" x 2.2" x 1.5"

Blade Length: 14.5"

Blade Diameter: 13.5Fr

Power Supply: 110/220 VAC; 50/60Hz

