

EVOLVE®

180W High Powered Diode



Case-by-Case
MOBILE SERVICE

The EVOLVE® Laser is used to vaporize and remove the obstructing prostate tissue that causes the symptoms associated with BPH. The 980nm wavelength provides the highest combined absorption in water and hemoglobin. Absorption of the laser energy by hemoglobin leads to excellent hemostasis. This 980nm laser allows for precise ablation of the obstructing tissue making it an ideal wavelength for this procedure.

Laser Specifications:

Laser type: Diode laser

Wavelength: 980nm

Power: up to 180W

Dimensions: 12" (h) x 22" (w) x 12" (d)

Weight: 66 lbs.



fortecmedical.com | 800.963.7101

✓ **Advanced Surgical
Technologies.**

✓ **Clinical Procedural
Support.**

✓ **Access for All.**