



**Case-by-Case**  
**MOBILE SERVICE**

# Single Use **Digital Ureteroscope**

.....

The Macrolux CoralView® flexible ureteroscope delivers high resolution imagery for detailed visualization and accurate navigation for effective treatment of stones.

The single use scope is individually packaged and pre-sterilized to eliminate scope reprocessing, maintenance, repair costs and any compromise in performance which are all common over the life of a reusable scope.

The scope has a reduced diameter of 7.5Fr providing improved irrigation and with a broad field of view of 120°. Its tip features a unique bullet-shaped profile for smoother insertion.



Reduce Thumb Fatigue: AutoHold technology  
self-locking lever with lightweight torque control

## [ THE FORTEC ADVANTAGE ]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

**Call today**  
**800.963.7101**  
to start ForTec Mobile Service

# Single Use Digital Ureteroscope

## Laser Specifications:

Deflection of distal tip: 275° up, 275° down

Diameter of insertion position: 7.5Fr


Field of View: 120°

Working length: 680mm

Instrumental channel: 3.6Fr


### Effortless access to small IPA calyces

275° deflection with  
instruments



### Enhanced stone-free rates

Easily drive flexible negative  
pressure sheath




### Excellent 1:1 torque rate

Soft and durable, resistant  
to damage



### Leak-proof

Patented T-valve with a  
water tap




### 7.5Fr. Superslim shaft by SlimEase Technology

Bullet-shaped tip for smoother insertion



### HD image quality by Xclear technology

Sharp image with accurate color reproduction



Full HD  
touchscreen



Video and voice  
recording



Anti-glare  
monitor



Fingerprint-free  
glass



FORTEC medical

>> Start **FORTEC** Mobile Service **Today**

[www.fortecmedical.com](http://www.fortecmedical.com) | 800.963.7101