



FORTECmedical

Digital Ureteroscope

Single Use



**Case-by-Case
MOBILE SERVICE**

The Macrolux CoralView® flexible ureteroscope delivers high resolution imagery for detailed visualization and accurate navigation for effective treatment of stones.

The single use scope is individually packaged and pre-sterilized to eliminate scope reprocessing, maintenance, repair costs and any compromise in performance which are all common over the life of a reusable scope.

The scope has a reduced diameter of 7.5Fr providing improved irrigation and with a broad field of view of 120°. Its tip features a unique bullet-shaped profile for smoother insertion.

Single-use form avoids the need for sterilization

A new scope every time ensures highest quality imaging

Slim diameter improves irrigation flow



Reduce Thumb Fatigue: AutoHold technology self-locking lever with lightweight torque control



**Advanced Surgical
Technologies.**



**Clinical Procedural
Support.**




Access for All.

Digital Ureteroscope Single Use


Effortless access to small IPA calyces

275° deflection with instruments



Enhanced stone-free rates

Easily drive flexible negative pressure sheath



Laser Specifications:

Deflection of distal tip: 275° up, 275° down

Diameter of insertion position: 7.5Fr

Field of View: 120°

Working length: 680mm

Instrumental channel: 3.6Fr


Excellent 1:1 torque rate

Soft and durable, resistant to damage



Leak-proof

Patented T-valve with a water tap



Full HD touchscreen



Video and voice recording



Anti-glare monitor



Fingerprint-free glass

7.5Fr. Superslim shaft by SlimEase Technology

Bullet-shaped tip for smoother insertion



HD image quality by Xclear technology

Sharp image with accurate color reproduction

