



FORTECmedical

PIRANHA

Morcellation System



Case-by-Case
MOBILE SERVICE

The PIRANHA Morcellation System offers a comprehensive, precisely coordinated system comprising a laser resectoscope, PIRANHA scope and morcellator for treatment of any size prostate or other morcellation procedure.

Powerful turbine suction rapidly draws in tissue for morcellation

Dynamic suction provides a clear field of view to maintain orientation during the case

Works with any brand of laser and fibers up to 1000 μ m

No need to remove sheath when changing instruments

Morcellation and suction controlled by one foot pedal



**Advanced Surgical
Technologies.**



**Clinical Procedural
Support.**



Access for All.

Product Features & Specifications

PIRANHA Morcellation System

Key Benefits

- High speed suction ensures maximum cutting force
- Single suction opening avoids inadvertent intake of the bladder mucosal membrane
- Oscillating cutting movement facilitates complete morcellation
- Morcellated tissue is big enough for histology
- Handle and cable are fully autoclavable



Operation

- Easy to toggle between different operating modes with the foot switch
- 'PIRANHA mode' for suctioning and morcellation
- 'Suction mode' for suctioning coagulated blood or very small tissue chips
- Tissue is cut, not chopped
- Disposable tube sets make work efficient

