



FORTEC_{medical}

UltraPulse CO₂ Laser System



**Case-by-Case
MOBILE SERVICE**

The UltraPulse CO₂ can be used for facial resurfacing for rhytids and acne scarring. The UltraPulse is also used for the removal of perioral and periorbital wrinkles, as well as blepharoplasty.

The UltraPulse CO₂ comes with the CPG (computer Pattern Generator). This helps treat large areas quickly with uniform fluence, allowing precise removal of tissue, layer by layer, with minimal thermal damage.



Laser Specifications:

Laser type: CO₂, Sealed cavity, RF excited, Tankless

Wavelength: 10.6 micron

Power: 1-100W

Dimensions: 77" (h) x 17" (w) x 33" (d)

fortecmedical.com | 800.963.7101



**Advanced Surgical
Technologies.**



**Clinical Procedural
Support.**



Access for All.