




FORTECmedical

Wolf Morcellator Tray Inventory

Images	Part Name	Part Number
	30° Telescope	8654.422
	Continuous Irrigation Outer Sheath (26Fr)	8655.374
	Continuous Irrigation Inner Sheath (24Fr)	8655.3841
	Visual Obturator	8415.12
	Working Insert/Laser Bridge	8632.284
	Working Insert with Continuous Irrigation	8654.911
	Stabilization Guide	8654.985
	Morcescope (12°)	8970402

Part Name	Rotatable Collar	Blue Bevel O-Ring	Yellow O-Rings (part of Morcescope)	Stopcock	Lauer Lock Connector
Part Number	8654.3742	15176.127	9500118	884.000	886.00
Images					
Part Name	Automatic Valve	Yellow Sealing Membrane	Blue Sealing Cap	Rubber Caps for Irrigation	Toomey Adapter (optional)
Part Number	8920.311	89.103	89.02	N/A	845.11
Images					
Part Name	Light Cord Adapter (30° scope)	Light Cord Adapter (to light source)	Light Cord Adapter (for Wolf Light Cord)	Light Cord Adapter (Morcescope)	Elick Adapter (optional)
Part Number	808801	8095.07	809509	8095.00	8858.14
Images					

Yellow bladder is disposable, but needs to be sterilized!

Part Name	Light Cord	Irrigation Adapter (part of Morcescope)	Laser Working Element
Part Number	80662530	15461.039	8654.282
Images			